

Manual Wheelchair Accessories



Wheelchair bags are available and can be placed on the back of the wheelchair or under the seat.

Cup or drink holders can be attached to carry drinks more easily.

A crutch or walking aid holder can be attached to some wheelchairs to carry the walking aids of more mobile individuals.



Oxygen bottles carriers can be hung over the back of the wheelchair. Care must be taken to ensure that the oxygen bottle does not unbalance the chair and cause it to tip backwards on slopes. An IV pole may also be attached to a wheelchair.



A range of arm supports can be placed directly on the armrest or used to replace the standard arm pad if the user requires more support for one arm.

Push mitts improve grip for self propulsion and will provide some protection to the hands.



A range of trays are available including flip up one sided trays and full width trays. These trays may be used to transport items or to assist in positioning the upper body.



Ankle and calf straps may be used to prevent the users foot from falling off the back of the footplate. A calf strap is particularly required if elevating leg rests are being used.

A range of belts and harnesses are available to assist in maintaining the users position in the chair, or for safety.

Rain capes are available to protect the user and wheelchair from the rain. Covers which go over the wheelchair are also available to protect the user from the weather.

NEED TO KNOW MORE?

Telephone, fax or e-mail your enquiry to obtain further information from one of the therapists about products. Alternatively, you can make an appointment to trial the equipment at our display centres or through our outreach service. Please phone LifeTec for further details.

Last Updated: March 2009